



Doctors Oxygen Service, Inc.
10180 South 54th Street
Franklin, WI 53132
800-558-3834
FAX: 414-421-7055

Apollo Anesthesia Machine - Test Equipment & Fittings

Part# 8602842 Fitting, NIST, O2



Part# 8602843 Fitting, NIST, N2O




Part# 8602844 Fitting, NIST, AIR



Part# 7910751 Fitting, M15 X 1




(3 each recommended – Screws into Above NIST Fittings)



Part# 7901297 Rubber Plug



Part# 7911268 Adapter, Pneumatic Interface



Part# M25647 Adapter, 22mm X 22mm



Part# 7911300 Port Puller Extraction Tool



Part# MP02400 Test Lung



Part# 1202537 Wacker Elastocil Silicone



Part# 1565559 Lubrication, Klüberalfa YV 93-302 4G



**Part# 7910349 Calibration Gas
(1% Sevo, 1% Iso, 5% CO₂, 70% N₂O, Bal O₂)**



Part# 7910346 Regulator, Calibration Gas



Part# TESTRIG-1 Leak Test Device, Tubing Asm.



Part# TESTRIG-2 Leak Test Device, Luer Lock Adapter
(Only needed If using Universal Pressure Meter W/ Luer Lock Connector)



Part# TESTRIG-3 Leak Test Device, Squeeze Bulb Adapter




Part# TESTRIG-5 Leak Test Device, Test Terminal Adapter
 (2 Each Required)



Part# TESTRIG-10 Leak Test Device, Gas Inlet Regulator Output
 Remove "T" fitting for calibrating Pressure Regulators PRPN2O, PRPAIR, PRPO2



Part# TESTRIG-11 Gas Manifold Test Device, Modified, Apollo



Part# M1-870-P Regulator, O2, "E" Cylinder Yoke
 (Adjusted for 70psi output) Used with Gas Manifold Above



Part# TESTRIG-12 Apollo Test Rig, Fresh Gas Measurement




Part# TESTRIG-13 **Apollo Test Rig, Regulator Performance Test**



Part# TESTRIG-14 **Apollo Test Rig, Downstream Pressure Test**



Part# TESTRIG-15 **Apollo Test Rig, PGM Gas Measurement**



Part# 2165619 **Breathing Hose 22mmX35cm (14 Inch)**



Part# 7911267 **Service Cable, Ventilator Extension Cable**
(For Diagnostic Testing of Ventilator when removed from Machine)



Part# 7901808 **Service Download Cable**
(Service Laptop to Apollo Anesthesia Machine Communication)



Part# 7901068 **SPO2 Shorting Cable**
(Only required if SPO2 Option is Installed)



Part# 7900894 **Hose Cutter**
(Used for proper cutting / trimming of Rigid Nylon Tubing)



Part# 7900953 **Socket 11mm – Modified**
 (Modified thin wall socket used for removal of Flow Control Valve Knobs)

Part# 7911785 **Flow Control Valve Removal Tool**



Part# 4070 **TSI Certifier W/High Flow Module & Carrying Case**
 For Accurate measurement of Gas flows from 15 – 300 L/min.



Part# 4072 **Low Flow Module, TSI Certifier**
 For Accurate Measurement of Gas Flows from 0.00 – 15.00 L/min.



Part# FI-8000P **Riken Gas Analyzer**
 (For Vaporizer Output Verification Testing)



Part# 417-138 **Fluke 179 Digital Multimeter W/Temperature Measurement**
 (Or Equivalent – Temp Measurement 0 - 100°C +/-1 °C)



Part# DPM-2001 **DPM-2000 Series Universal Pressure Meter**
 (Or Equivalent – Measurements up to 75psi)



Part# 420-659 Dale 601 Electrical Safety Analyzer
(Or Equivalent)



Service Laptop Computer Obtain Locally

Used for accessing Digital Service Manuals & Documentation (Draeger ServiceConnect 2.0)
Running Maintenance and Testing Software Applications (Draeger Service Software)
Downloading of service logs
Loading and/or Updating Machine Software