




Doctors Oxygen Service, Inc.
9760 South 60th Street
Franklin, WI 53132
800-558-3834
FAX: 414-421-7055

Fabius / TIRO Anesthesia Machine - Test Equipment & Fittings

Part# TESTRIG-1 Leak Test Device, Tubing Asm.



Part# TESTRIG-2 Leak Test Device, Luer Lock Adapter



Part# TESTRIG-3 Leak Test Device, Squeeze Bulb Adapter



Part# TESTRIG-4 Leak Test Device, Fabius Fresh Gas Adapter



Part# TESTRIG-5 Leak Test Device, Test Terminal Adapter
2 each Required



Part# TESTRIG-6 Leak Test Device, 2 Foot Extension

(2 Each)



Part# TESTRIG-7 Leak Test Device, Pipeline Check Valve Leak




Part# TESTRIG-8 Leak Test Device, Control Port Leakage



Part# TESTRIG-9 Leak Test Device, Airway Pressure & APL Port Connections




Part# TESTRIG-10 Leak Test Device, Gas Inlet Regulator Output



Part# 2165619 Hose, 22mm X 12 inch



Part# TESTRIG-11 Gas Manifold Test Device, Modified, Apollo



Part# M1-870-P Regulator, O2, "E" Cylinder Yoke
(Adjusted for 70psi output) Used with Gas Manifold Above



Part# MP02400 Test Lung



Part# 1202537 Wacker Elastocil Silicone



Part# 1565559 Lubrication, Klüberalfa YV 93-302 4G



Part# 7900894 Hose Cutter
(Used for proper cutting / trimming of Rigid Nylon Tubing)



Part# 7911785 Flow Control Valve Removal Tool



Part# 7900953 Socket 11mm – Modified
(Modified thin wall socket used for removal of Flow Control Valve Knobs)

Part# 4070 TSI Certifier W/High Flow Module & Padded Carrying Case
For Accurate measurement of Gas flows from 15 – 300 L/min.



Part# 4072 Low Flow Module, TSI Certifier
For Accurate Measurement of Gas Flows from 0.00 – 15.00 L/min.



Part# FI-8000P **Riken Gas Analyzer**
 (Or Equivalent – For Vaporizer Output Verification)



Part# 417-138 **Fluke 179 Digital Multimeter**
 (Or Equivalent)



Part# DPM-2001 **DPM-2000 Series Universal Pressure Meter**
 (Or Equivalent)



Part# 420-659 **Dale 601 Electrical Safety Analyzer**
 (Or Equivalent)



Service Laptop Computer **Obtain Locally**
 Used for accessing Digital Service Manuals & Documentation (Draeger ServiceConnect 2.0)
 Running Maintenance and Testing Software Applications (Draeger Service Software)
 Downloading of service logs
 Loading and/or Updating Machine Software