


AGR

ANESTHETIC GAS RECLAMATION

21st Century Scavenging




T
O
D
A
Y




NOTE: AGSS VACUUMS CONTINUOUSLY DRAWS REGARDLESS OF ON/OFF POSITION

T
O
D
A
Y



BENEFITS OF DGSS

- Vacuum Pump Runtime Cut in Half
- Typically 1-2 year investment payback
- Can be used on future anesthesia machines (i.e. do not have to purchase a scavenger with a new machine)
- Sets up for future gas reclamation



VAC PUMP HOUR METER

0 0 7 4 3 2





:

0 0 7 6 4 7

NOW TYPICALLY 7-8 HOURS
PER PUMP PER DAY

T
O
M
O
R
R
O
W

ANESTHETIC MACHINE



NOW TYPICALLY 7-8 HOURS PER PUMP PER DAY

NOTE: DGSS ONLY DRAWS VACUUM WHEN PATIENT ON TABLE AND PRESSURE SWITCH TRIPPED